



Herbert Wertheim
College of Engineering
UNIVERSITY of FLORIDA

CAMbox


Compact Autonomous Microbioreactor



Section 22695, Group 14

Marcela Abadia, Joao Pedro Dos Santos, Ashton Goan, Ollie Goodall, Brenden Modi, Audrey Myat, Wojciech Przepiora, Pelayo Urrios



- **Compact:** Smaller footprint allows for easier placement.
- **Lower cost:** Fewer materials used enabled use of higher quality OTS parts.
- **Marketable:** Compact, cost efficient, design makes it available to a larger audience.





Enclosure

M3A: Each unassembled component weighs under 50 pounds

C3: Moveable by one person after disassembly

M4A: Less than 36 inches wide

M4B: load of assembly is lower than 600 pounds

M4C: Bottom face dimensions must be smaller than 24 inches wide and 60 inches long

C4: Fits on a research benchtop,

M6: Access hole must allow two hands to fit through at least 20 inches long, 5 inches wide

C6: Has an easily accessible interior for cleaning

M12: Signs will be posted to the device warning users about volatile materials and to keep them far from the reactor

C12: Only nonreactive materials contact lab chemicals

M13: Waste will be stored and made safe until waste can be fully disposed of


C13: Appropriate for operation in a BSL-2 space

M14: Exterior casing must always remain at a temperature under 44 Celsius or 111 Fahrenheit

C14: Has an exterior surface that is not too hot to comfortably touch




Enclosure

- **Material Selection**
 - Constructed from ASTM A36 Steel
 - Layer of Polyurethane Foam Rigid
 - Insulating R-value ≈ 5
 - Maintains internal temperatures



Enclosure

- Material Analysis
 - ASTM A36 Steel has a yield strength of $2.5e8$ Pa
 - Maximum stress stress from gantry on supports $\approx 1.271e5$ Pa



Liquid Handling

M11: Disposable tips will be made of nonporous materials. 1 disposable tip per liquid being used.

C11: Only nonporous materials contact cell cultures

M15: Waste will be neutralized by being placed in a reservoir containing 0.10 of the volume in bleach.

C15: Is capable of sequestering and neutralizing its own liquid and solid waste

M34: Quantified by the number of steps the subsystems require to complete the process of adding or subtracting liquid. The number of steps this subsystem takes will be between 3 and 6.

C34: Is capable of automated liquid handling with fluid addition/subtraction from each well or tube

M35: To determine whether aerosols will be created the Reynolds number will be calculated. The Reynolds number will have a maximum of 2300 in order to avoid creating aerosols

C35: Dispenses fluid without creating aerosols

M36: Dispenses fluid from 225 $\mu\text{L/s}$ to 300 $\mu\text{L/s}$

C36: Achieves dispense rates from 225 $\mu\text{L/s}$ to 300 $\mu\text{L/s}$

M37: Deposit from 5-20000 μL of aliquot fluid volume

C37: Deposits a minimum/maximum aliquot fluid volume from 5-20,000 μL

M38: Volume accuracy 0.01 μL

C38: Achieves dispensing volume accuracy of $\pm 0.1 \mu\text{L}$




M39: Volume accuracy 0.01 μL

C39: Achieves dispensing volume precision $\pm 0.01 \mu\text{L}$

M40: Cross contamination will be avoided by using disposable tips. The number of tips used will depend for the number of pipettes, a minimum of 1.

C40: No cross contamination between individual wells/tubes during liquid handling

Liquid Handling



Liquid handling

- Flow rate: 375 $\mu\text{L}/\text{s}$

$$\bar{V} = \frac{2\pi}{60} \times 2250 \times (9.566 \times 10^{-3}) = 2.25 \text{ m/s}$$

$$Q = 2.25 \times (1.66 \times 10^{-7}) = \frac{3.75 \times 10^{-7} \text{ m}^3}{\text{s}} = 375 \mu\text{L}/\text{s}$$

- No Aerosols, confirmed with Reynolds Number:

$$Re = \frac{997 \times 2.25 \times 4.6 \times 10^{-4}}{1.08 \times 10^{-3}} = 957$$

- Dispense accuracy: $\pm 0.03 \mu\text{L}$

$\pm 0.05^\circ \rightarrow 7200 \text{ steps per revolution}$

screw lead: 1.5 mm

$$\frac{1.5}{7200} = 0.000208 \text{ mm/step}$$

Cross sectional area of syringe: $\frac{\pi(9.57 \times 10^{-3})^2}{2} = 0.00014 \text{ m}^2$

$$0.000208 \times 0.00014 = \pm 3 \times 10^{-11} \text{ m}^3$$

$$\pm 0.03 \mu\text{L}$$

User Interface

M7: Quantified by how long it takes to stop all function, a minimum of 0.5 seconds will be used and a maximum of 1 second

C7: Includes an easily actuated emergency shut-off that safely stops all functions

M8: There will be a minimum of 1 setting per number of subassemblies of the system and a maximum of 2 settings per number of subassemblies, making the minimum settings on the interface 8 and the maximum 16

C8: Has an intuitive user interface

M9: This need will be quantified by the brightness of the interface display or flashing of the buttons on a scale of 1-5, 1 being the dimmest and 5 being the brightest



C9: Has a visual indicator that is easily seen by the user and nearby personnel

M10: Each subsystem can be controlled individually, quantified by the number of controls on the interface, a minimum of 8 to allow control of different subsystems

C10: Is programmable: control parameters can be changed by the user, and more complex processes can be added

User interface

- Arduino will be inside the control box.
- Device connects to laboratory computer
- Sensors are compatible with Arduino.



Climate Control

M5: Operate with a 120V AC

C5: Runs from a single standard 120 VAC wall outlet

M17: Independently control at least 3 separate environments

C17: Maintains environmental conditions independently for each well plate or tube

M20: 4-70 degrees Celsius

C20: Maintains cultures in a well plate or tube at a constant temperature within a range from 4C to 70C

M21: 2.5 degree difference

C21: Maintains internal setpoint temperature with time and spatial variation less than $\pm 2.5C$

M22: 5 degree difference

C22: Uniformly heats/cool wells/tubes within the desired temperature range

M23: Covers less than .10 of surfaces on well plates or conical tubes


C23: Mitigates condensation on cooled well plate and tube surfaces

M24: 15 minutes

C24: Reaches setpoint culture temperature in less than 15 minutes


Climate Control - Cooling

- **Seifert 3152303 thermoelectric cooler**
 - More compact
 - Better longevity
 - No moving parts
 - 200W cooling capacity
 - Temperatures as low as -10 C



Climate Control - Cooling

- Compact enclosure enables rapid heating & cooling.
- Variable-speed fan offers greater user control.
- 20C to 4C in under 12 minutes.



Climate Control – Circulation

- System utilizes forced convection.
- Variable-speed Supermicro fan.
- 25.2 CFM with an outlet of 1.57in x 1.57in.



Maximum output air speed calculations:

Outlet area:

$$1.57[\text{in}] * 1.57[\text{in}] * \left(\frac{2.54[\text{cm}]}{1[\text{in}]}\right)^2 * \left(\frac{1[\text{m}]}{100[\text{cm}]}\right)^2 = 1.5903 * 10^{-3}[\text{m}^2]$$

Volumetric flow rate conversion:


$$25.2 \left[\frac{\text{ft}^3}{\text{min}}\right] * \left(\frac{1}{3.28}\right)^3 \left[\frac{\text{m}^3}{\text{ft}^3}\right] * \frac{1}{60} \left[\frac{\text{min}}{\text{s}}\right] = 0.0119 \left[\frac{\text{m}^3}{\text{s}}\right]$$

Maximum output air speed:

$$\frac{0.0119 \left[\frac{\text{m}^3}{\text{s}}\right]}{1.5903 * 10^{-3}[\text{m}^2]} = 7.48 \left[\frac{\text{m}}{\text{s}}\right]$$


Climate Control - Heating

- 200W DBK FGC3000 Series Fan Heater.
- Built in DIN rail clip.
- Maximum temperature of 70 C.
- 300,000h+ lifetime.
- Secondary function as defogger during cooling




Climate Control - Heating

Air Velocity v. Time to Reach Setpoint –
Heating of a 50mL Conical Tube




- With selected fan: 20C to 70C in under 14 minutes.

- Analysis based on forced convection: cylinder in cross-flow



Climate Control - Heating



- Areas, radii, lengths from conical tube dimensions:

- $A_1 = \pi d_o h = \pi * 30\text{ [mm]} * 115\text{ [mm]} = 10,838.49\text{ [mm}^2\text]}$
- $A_3 = \pi d_i h = \pi * 26\text{ [mm]} * 115\text{ [mm]} = 9393.36\text{ [mm}^2\text]}$
- $r_1 = 13\text{ [mm]}, r_2 = 15\text{ [mm]}, L_2 = 115\text{ [mm]}, L_3 = 13\text{ [mm]}$

- Thermo properties for polypropylene (conical tube) surface at 20 C, air at 70 C, water at 20 C:

- $v_{air} = 20.92 * 10^{-6}\left[\frac{\text{m}^2}{\text{s}}\right], k_{water} = 598.03\left[\frac{\text{W}}{\text{m}\cdot\text{K}}\right], k_{air} = 30 * 10^{-3}\left[\frac{\text{W}}{\text{m}\cdot\text{K}}\right]$
- $Pr = 0.7, Pr_s = 0.707, k_{PP} = 0.165\left[\frac{\text{W}}{\text{m}\cdot\text{K}}\right]$
- For $1000 < Re < 2 * 10^5$: $C = 0.26, m = 0.6$
- For $Pr < 10$: $n = 0.37$

- Modified convective heat transfer coefficient:

- $\bar{h} = \frac{\overline{Nu}_d k}{D}$

- $\overline{Nu}_d = C Re_D^m Pr^n (Pr \setminus Pr_s)^{\frac{1}{4}}$

- $Re_D = \frac{vD}{\nu}$

- For fan speed 7.48 m/s:

- $Re_D = \frac{7.48\left[\frac{\text{m}}{\text{s}}\right] * 0.03\text{ [m]}}{20.92 * 10^{-6}\left[\frac{\text{m}^2}{\text{s}}\right]} = 10,727$

- $\overline{Nu}_d = 0.26(10727)^{0.6} (0.700)^{0.37} (0.700 \setminus 0.707)^{\frac{1}{4}} = 59.695$

- $\bar{h} = 59.695 * \frac{30 * 10^{-3}\left[\frac{\text{W}}{\text{m}\cdot\text{K}}\right]}{0.03\text{ [m]}} = 59.695\left[\frac{\text{W}}{\text{m}^2\cdot\text{K}}\right]$

- Heat transfer rate:

- $q = \frac{\Delta T}{R_{tot}} = \frac{343\text{ [K]} - 293\text{ [K]}}{2.748\left[\frac{\text{K}}{\text{W}}\right]} = 18.197\text{ [W]}$

- Time required to reach setpoint:

- $t = \frac{Q}{q} = \frac{14974\text{ [J]}}{18.197\text{ [J/s]}} = 789.31\text{ [s]} = 13.16\text{ [min]}$

Gas Control

M1: 10 years

C1: Operational lifetime of at least 10 years


M18: Inject 5 gases

C18: Is capable of safely injecting, measuring, and regulating the composition of the following gases in each compartment holding a well plate or tube: N₂, O₂, CO₂, CH₄, H₂



M19: 2 weeks

C19: Is capable of incubation periods long enough to permit 1,000 E. coli culture generations (~two weeks)

Gas Control



- Solenoid valves ASCO Red Hat
- Long life Equipment (5-20 million cycles)
- Two weeks of incubation





A mind map diagram with a central node labeled "OD/FI" in a rounded rectangle. Six lines radiate from this central node to six pairs of text. Each pair consists of a requirement (M) and a constraint (C). The lines are color-coded: purple for M2, M27, M30, and M31; green for M32; and blue for M33. The constraints (C) are connected to their respective requirements by horizontal lines.

OD/FI

M2: \$10,000

C2: Prototype production cost does not exceed \$10,000

M27: Illuminate from 400nm-640nm at 1kW/cm²

C27: Includes a photobioreactor mode to illuminate photosynthesis-capable cell cultures (e.g., cyanobacteria) with white light

M30: Measure OD at 600nm

C30: Measures optical density (OD) in all individual wells and conical tubes

M31: Excite fluorescence up to 640nm and detect emission up to 740 nm

C31: Measures fluorescent intensity (FI) in all individual wells and conical tubes


M32: Maintain Light intensity of 640nm

C32: . For OD/FI measurement, sustains adequate light intensity to make measurements at wavelengths not lethal to cells

M33: 6.5 Minutes


C33: Processes a 384 well plate through OD/FI measurements in less than 6.5 minutes

OD/FI Reader



- Includes:
 - PMT Detector
 - Filter Cube Chassis
 - Excitation Filter glass
 - Emission Filter Glass
 - Dichroic Mirror
 - Standard Mirror
 - Xenon Lamp

- ▶ Absorbance (OD)
- ▶ FI Excitation
- ▶ FI Emission
- ▶ Xenon Excitation



■ Analysis Speed

$$\frac{\text{Well Plate Dia.} * (\text{columns} * \text{rows} + 2)}{\text{Gantry Velocity Speed} * \text{Load Percentage}} + (\text{Analysis} * \text{well plates})$$

$$= \frac{0.0031 * (24 * 16 + 2)}{.513 * .6} + (0.15 * 384) = \mathbf{61.95 \text{ Seconds}}$$

■ Cost

- OTS Parts: \$1507.02
- Material: \$53.68
- Manufacturing: \$60.30

OD Capability

- Measures absorbance at 600 nm

FI Capability

- Measures Texas Red Fluorescence

White Light

- Illuminates well plate sample at 640nm 0.2 Kw/cm2

Lethality

- The 600nm featured ranges ensure all cells will live!

Well-Plate Shaker

M16: Every compartment is interchangeable — C16: Cultures microbes in fully enclosed compartments or vessels that are interchangeable

M29: Adaptable to shake well plates and tubes — C29: Has shaking patterns that are independent for each well plate or tube

M28: Capable of three shaking patterns — C28: Shakes well plates and tubes in Linear, Orbital, and Double Orbital patterns

M25: Hold well plates of 6,24,48,96, deep 96, 385 — C25: Accommodates existing culture well plates of the following sizes: 6, 24, 48, 96, deep 96, 384

M26: Hold conical tubes of 15mL and 50mL — C26: Accommodates existing conical tubes of the following sizes: 15mL & 50mL

Shaker

- Linear, orbital, and double orbital shaking patterns
 - Uses springs and two motors to achieve a maximum velocity of 14.4963 m/s
- Interchangeable between well plates and tubes
- Tray designed with non-porous material


Governing Equations



$$\omega = \frac{RPM}{60 \text{ s/min}} * 2\pi \frac{rad}{rev}$$



$$\omega = \frac{600 \text{ RPM}}{60 \text{ s/min}} = 62.83 \frac{rad}{s}$$

$$v = \omega r = 62.83 \frac{rad}{s} * 0.0381 \text{ m} = 2.39 \frac{m}{s}$$


$$v = \frac{10900 \text{ RPM}}{60} * 2\pi * 0.012 \text{ m} = 14.4963 \frac{m}{s}$$







Shaker – Acceleration



Acceleration to cause spillage and cross contamination of cell cultures

$$\frac{\partial P}{\partial x} = -\rho \cdot a_x = 0 \rightarrow \text{No acceleration in the } x \text{ axis}$$

$$\frac{\partial P}{\partial y} = -\rho \cdot a_y$$

$$\frac{\partial P}{\partial z} = -\rho \cdot (a_z + g) = -\rho \cdot g \rightarrow \text{No acceleration in the } z \text{ axis}$$

$$dP = \frac{\partial P}{\partial x} dx + \frac{\partial P}{\partial y} dy + \frac{\partial P}{\partial z} dz = -\rho \cdot a_y \cdot dy - \rho \cdot g \cdot dz = -\rho \cdot (a_y \cdot dy + g \cdot dz)$$



$$\frac{\partial P}{\partial z} = -\rho \cdot g \rightarrow \int_{P_2}^{P_1} \partial P = \int_{Z_2}^{Z_1} -\rho \cdot g \cdot dz \rightarrow P_1 - P_2 = -\rho \cdot g(Z_1 - Z_2)$$

Apply equation along surface of the liquid, since pressure is constant

$$dP = 0 \rightarrow \text{Constant pressure in the surface}$$

$$0 = -\rho \cdot (a_y \cdot dy + g \cdot dz)$$

$$a_y \cdot dy = -g \cdot dz \rightarrow -\frac{a_y}{g} = \frac{dz}{dy}$$



Shaker – Acceleration (Cont.)

Relationship of y since volume is constant

$$y = 2 \frac{h_o L}{H}$$




$$a_y = -g \cdot \frac{Z_2 - Z_1}{Y_2 - Y_1} = -g \cdot \frac{z}{y} = -g \cdot \frac{z \cdot H}{2h_o L}$$

$$P_2 = \rho g (Z_1 - Z_2)$$

Thus, maximum acceleration before the cell cultures overspill is

$$y = 2 \frac{h_o L}{H} = 2 \frac{(4 \text{ mm})(6.49 \text{ mm})}{10.3 \text{ mm}} = 5.39 \text{ mm}$$

$$a_y = -g \cdot \frac{z}{y} = -9.81 \frac{\text{m}}{\text{s}^2} \left(\frac{10.3 \text{ mm}}{5.39 \text{ mm}} \right) = 18.745 \text{ m/s}^2$$



Shaker – Acceleration (Cont.)

Maximum acceleration of the shaker:

Linear:

$$v = \frac{10900 \text{ RPM}}{60} * 2\pi * 0.012 \text{ m} = 14.4963 \frac{\text{m}}{\text{s}}$$
$$\alpha = \frac{\pi(10900 \text{ RPM})}{1(30)} = 1141.445 \frac{\text{rad}}{\text{s}^2}$$
$$a = \alpha r = 1141.445 \frac{\text{rad}}{\text{s}^2} * 0.012 \text{ m} = 13.69 \frac{\text{m}}{\text{s}^2}$$

Orbital and Double Orbital:

$$\alpha = \frac{\pi(13900 \text{ RPM})}{1(30)} = 1455.6$$
$$a = \alpha r = 1455.6 \frac{\text{rad}}{\text{s}^2} * 0.012 \text{ m} = 17.47 \frac{\text{m}}{\text{s}^2}$$

Cost

Subsystem	Total cost
Enclosure	\$1468.54
Liquid handling	\$323.47
User interface	\$243.99
Climate control	\$1903.44
Gas control	\$1071.85
OD/FI	\$1580.61
Shakers	\$343.20
Total	\$6074.39

Why CAMbox should be selected for prototyping

- Efficient use of space and funds.
- Possibility to market outside of universities.
- All customer needs were met.

Conclusion

- CAMbox is a compact, cost efficient design.
- The autonomous microbioreactor is a single environment and can house cell cultures for 2 weeks.
- This design is meant to be targetted at an audience outside of the University.

Thank you!



UF | Herbert Wertheim
College of Engineering
UNIVERSITY *of* FLORIDA

POWERING THE NEW ENGINEER TO TRANSFORM THE FUTURE